

\\INVESTIGATORS U-Z\IVANELZAKKER\CCJ pilot\protocol setup\ContinuousStim_ep2d_bold_1p5m m_iso_R3_SMS4_5x5

TA: 7:15 PM: REF Voxel size: 1.5×1.5×1.5 mmPAT: 12 Rel. SNR: 1.00 : epbMp12

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Resolution - iPAT

Accel. factor PE	3
Ref. lines PE	48
Accel. factor slice	4
Reference scan mode	EPI/separate

Resolution - Filter Image

Distortion Corr.	Off
------------------	-----

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
Hamming	Off

Routine

Slice group	1
Slices	92
Dist. factor	0 %
Position	R1.1 A16.6 H0.0 mm
Orientation	T > S2.7 > C0.8
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	1190 ms
TE	22.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	701-732;H12-18,21-28,31-38,41-44

Geometry - Common

Slice group	1
Slices	92
Dist. factor	0 %
Position	R1.1 A16.6 H0.0 mm
Orientation	T > S2.7 > C0.8
Phase enc. dir.	A >> P
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	1190 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	R1.1 A16.6 H0.0 mm
Orientation	T > S2.7 > C0.8
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A1.0 H1.8
Phase	-1.0 mm
Read	0.0 mm
Shift	1.8 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Contrast - Common

TR	1190 ms
TE	22.0 ms
MTC	Off
Flip angle exc	60 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	326
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

Accel. mode	Slice accel.
-------------	--------------

Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Performance
AutoAlign	Head > Brain
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	R1.1 A16.6 H0.0 mm
Orientation	T > S2.7 > C0.8
Rotation	1.97 deg
A >> P	192 mm
R >> L	192 mm
F >> H	138 mm
Reset	Off

System - Tx/Rx

Frequency 1H	297.189966 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	1190 ms
Concatenations	1

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Motion correction	Off
Spatial filter	Off
Measurements	326
Delay in TR	0 ms
Multiple series	Off

Sequence - Part 1

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.64 ms

Sequence - Part 1

Bandwidth	1860 Hz/Px
-----------	------------

Sequence - Part 2

EPI factor	128
RF pulse type	Low SAR
Gradient mode	Fast

Sequence - Special

Imaging Dummy TRs	3
ACS Dummy TRs	3
FLEET Dummy Pulses	5
SMS Shift	3
RF Clip	0
VERSE Factor	1.75
ACS Type	ALL FLEET
FLEET FA	10 deg
Physio recording	Off
Wait for TCP/IP Trigger	Off
Reverse Phase Encoding	Off

\\INVESTIGATORS U-Z\IVANELZAKKER\CCJ pilot\protocol setup\BlockStim_ep2d_bold_1p5mm_iso_R3_SMS4_5x5

TA: 7:15 PM: REF Voxel size: 1.5×1.5×1.5 mmPAT: 12 Rel. SNR: 1.00 : epbMp12

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Resolution - iPAT

Accel. factor PE	3
Ref. lines PE	48
Accel. factor slice	4
Reference scan mode	EPI/separate

Resolution - Filter Image

Distortion Corr.	Off
------------------	-----

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
Hamming	Off

Routine

Slice group	1
Slices	92
Dist. factor	0 %
Position	R1.1 A16.6 H0.0 mm
Orientation	T > S2.7 > C0.8
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	1190 ms
TE	22.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	701-732;H12-18,21-28,31-38,41-44

Geometry - Common

Slice group	1
Slices	92
Dist. factor	0 %
Position	R1.1 A16.6 H0.0 mm
Orientation	T > S2.7 > C0.8
Phase enc. dir.	A >> P
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	1190 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	R1.1 A16.6 H0.0 mm
Orientation	T > S2.7 > C0.8
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A1.0 H1.8
Phase	-1.0 mm
Read	0.0 mm
Shift	1.8 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Contrast - Common

TR	1190 ms
TE	22.0 ms
MTC	Off
Flip angle exc	60 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	326
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

Accel. mode	Slice accel.
-------------	--------------

Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Performance
AutoAlign	Head > Brain
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	R1.1 A16.6 H0.0 mm
Orientation	T > S2.7 > C0.8
Rotation	1.97 deg
A >> P	192 mm
R >> L	192 mm
F >> H	138 mm
Reset	Off

System - Tx/Rx

Frequency 1H	297.189966 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	1190 ms
Concatenations	1

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Motion correction	Off
Spatial filter	Off
Measurements	326
Delay in TR	0 ms
Multiple series	Off

Sequence - Part 1

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.64 ms

Sequence - Part 1

Bandwidth	1860 Hz/Px
-----------	------------

Sequence - Part 2

EPI factor	128
RF pulse type	Low SAR
Gradient mode	Fast

Sequence - Special

Imaging Dummy TRs	3
ACS Dummy TRs	3
FLEET Dummy Pulses	5
SMS Shift	3
RF Clip	0
VERSE Factor	1.75
ACS Type	ALL FLEET
FLEET FA	10 deg
Physio recording	Off
Wait for TCP/IP Trigger	Off
Reverse Phase Encoding	Off